



## Medical Filtration Solutions Ltd.

Company Registration Number: 8718303


72 Roman Way Industrial Estate,  
Longridge Road,  
Preston,  
Lancashire.  
PR2 5BE  
Tel: 01772 790820  
Fax: 01772 797515

### TECHNICAL DATA SHEET

**Product Code - M01.4.008**

**Product Description - S8 Refill**

Media	Polyester fibre . class G4
Compatibility	Compatible with ResMed S8 Compact, Elite, Autoset Vantage & Escape CPAP devices.
Alternatives	Other Efficiency Variants also Available



DIMENSIONS ARE APPROXIMATE

MFS-TDS- 0015  
Date: 24-12-2014

QF: 7.3-18  
Rev: A