

SPIRO *flo*



Medical Filtration Solutions Ltd.

Company Registration Number: 8718303



72 Roman Way Industrial Estate,
Longridge Road,
Preston,
Lancashire.
PR2 5BE
Tel: 01772 790820
Fax: 01772 797515

TECHNICAL DATA SHEET

Product Code - M02.2.015

Product Description - Spirometry filter & Bitegrip

Filter media	Statically charged filtration media
Bacterial Filtration Efficiency	99.9999%
Viral Filtration Efficiency	99.999%
Resistance to flow @ 720 L/min	63 Pa
Dead Space	71 ml
Material	Polypropylene
Connections	Comfortable oval mouthpiece for patient use. Device fitment to suit: Jaeger - Masterlab, Masterscreen & Masterscope series. Ferraris Medical - PK Morgan - Autolink-Benchmark, Transferflo and USA Transfer systems. Sensormedics - 2200-2400-2450 and Vmax systems. & Pneumotach flowmeters. Note: Some models may require adapters.



Dimensions are approximate

Date: 04-02-2016
Issue: 1

QF: 7.3-18
Iss: 4